
ClO_2 GENERATOR

INSTRUCCIONES DE USO
USER MANUAL
GEBRAUCHSANWEISUNG
NÁVOD K POUŽITÍ
HASZNÁLATI ÚTMUTATÓ
INSTRUKCJA UŻYCIA
MANUEL DE L'UTILISATEUR

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Distribuidor:

MEMBRANE ELECTROLYTIC GENERATOR OF ULTRA-PURE CHLORINE DIOXIDE

Thanks for purchasing a chlorine dioxide generator. This kind of generator produces an ultrapure gas called chlorine dioxide (ClO_2) using a special type of electrolysis and a two-phase semi-permeable membrane.

The ClO_2 generator has been designed for laboratory use. Attention! Chlorine dioxide is an explosive and flammable gas! As a matter of principle, only carry out production in outdoor areas with sufficient air flow, away from fire or spark sources, etc. Avoid exposing the generator to sunlight!

Strictly observe the safety instructions during manufacture.

The manufacturer is not responsible for any damages caused by improper use of the generator.

The device only is designed for the production of ClO_2 from sodium chlorite solution in the concentration range of 20-40%. Never use any other solution!

The high-quality food grade steel, pure platinum and other durable materials have been used for the production of electrodes. Nevertheless, a slight corrosion of the electrode surface may occur, as the sodium chlorite is strongly active due to its pH. However, this phenomenon has no impact on the purity of the final product - ultrapure chlorine dioxide. The purity of the ClO_2 produced by this generator exceeds 99.99%. It is practically impossible to achieve such a purity of ClO_2 by chemical means.

ClO_2 purity

ClO_2 from this generator is up to 500 times purer than ClO_2 produced chemically. The purity of the produced ClO_2 can be indirectly measured with a conductivity meter. Detailed instructions for this measurement have been posted on www.clo2gen.com for the inspection of any ClO_2 generator user.

You now have a completely unique instrument, which has been developed over a period of five years by several top European specialists.

Safety instructions:

Always read these safety instructions before each use of the generator!

- Any improper use of this product may result in its damaging!
- Never use the generator near a source of open flame!
- Place and keep the generator completely out of the reach of children. Place and keep the generator out of reach of persons who have not been instructed regarding the generator, its use and its function!

- Wear protective equipment!
- Wear safety glasses!
- Wear rubber protective gloves!
- Never inhale ClO_2

- Do not eat or drink while the generator is running!

- The generator is delivered without liquid. However, after the first use, the generator must be continuously filled with water (during storage) or sodium chlorite (24.5% solution) during production. This is necessary in order to maintain the tightness of the semi-permeable membrane inside the generator. The generator can be drained for a short period of time (12-15 hours), but it is necessary to close both filling holes immediately after draining and then seal the generator with the included silicone tubing - hereafter referred to as „travel condition“. This will maintain the 100% humidity inside both chambers of the generator. Failure to observe this essential condition may result in irreversible damage to the generator.

Travel condition

- Protect the generator from frost!
- protect the generator from direct sunlight (especially when making ClO_2 !)
- if the unit is out of service for a long time, store it in a dark place and filled with water.
- Never switch on the ClO_2 generator at the same time as a UV source or ozone generator. Do not move the generator once it has been put into operation. **Never mix the solutions in the two chambers** In an event the generator still needs to be transported, drain the generator cartridges according to the instructions and fill the generator with clean water for transport.

Preparation for first chlorine dioxide production:

Place the generator in a horizontal location in an open or well-ventilated place.

Follow the images and text below.

Travel condition

Unscrewed, ready.

Insert the funnel into the first chamber and pour sodium chlorite up to the line on the generator.

Poured up to the mark. Good.


Then fill the second chamber of the generator with sodium chlorite.


Prepare the generator, funnel and sodium chlorite solution.


Unscrew the two holes to pour the sodium chlorite into the generator chambers.


Screw on the caps as shown and tighten securely.


Do not replace the caps! The outlet marked ClO_2 serves for the insertion of the tubing!


Unscrew the connector cover for the electrode connection.


Connect the connector to the generator and screw in the lock.


Screwing the connector lock.


Connect the electronics module to the generator cable.


Consequently, connect the air compressor to the generator - to the outlet labelled AIR hose with the one-way valve. Afterwards connect the little air compressor to the USB socket on the electronics module.

Connect the power supply - mains adapter or power bank.


Attach a hard tube to the outlet labelled CIO2 and insert the other end into the water bottle, all the way to the bottom.

Potentiometer - Setting the amount of bubbles to carbonate the solution

Everything is ready, now. Switch on the power and ClO_2 production will start in a few minutes.

Use the potentiometer on the electronics module to set the optimum carbonation of the water with ClO_2 bubbles. In the beginning of the production, slower "bubbling" is advisable; after about 1 hour move to a faster bubbling. Bubbles must ALWAYS come out of the generator when it is on. Otherwise, even though you have tried different potentiometer positions, turn the generator off, check the wiring and if ClO_2 bubbles are not discharged into the water bottle, turn the generator off and contact your dealer.

Start of production. The ClO_2 chamber begins to colour.

Generator filling colour at full power

We recommend using chilled distilled water for carbonation. Full power and production of ultrapure chlorine dioxide with a new generator cartridge is reached after about 20 minutes of operation. Until then, ultrapure chlorine dioxide is produced with less intensity, as impurities in the sodium chlorite are eliminated. After the sodium chlorite has been purified the generator runs at full power.

Use up the generator charge at once. Make a supply of ClO_2 solution for a longer time Store the produced solution in a cool and dark place.

If you follow these procedures, the generator will serve you flawlessly for many years.

When the solution capacity is used up (the sodium chlorite solution „fades“ during the operation), disconnect the generator from the power supply. Disconnect the electrode power cord from the generator, screw in the connector plug, disconnect the tubing, and pour the fillings of both chambers into a pre-filled bucket of water. Pour the contents of the bucket down the drain. Then rinse both chambers of the generator at least twice with clean water.

Then fill the generator with water, close both chambers and connect both outlets with a hose to ensure a good seal and store out of reach of children. Leave the generator like this until a further production.

The complete exhaustion of the generator's filling capacity is indicated by the fact that the new clean water ready for ClO_2 saturation will no longer be yellow and the solution inside the generator will have significantly lighter colour.

The inner membrane of the generator must be kept moist at all times, otherwise it will deteriorate the generator!

Only use outdoors!

Never use near fire!

Wear protective equipment including goggles!

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